

LIPOSORBER®

Provides Hope When Drug Therapy Fails™

For Your Pediatric and Adult Primary FSGS Patients

Enabled nearly **50%** of steroid-resistant primary FSGS patients to attain complete or partial remission, 2 years after treatment.²

Consider LIPOSORBER therapy for primary FSGS patients when:


- standard therapies have failed or are not well-tolerated & patient's $GFR \geq 60 \text{ ml/min/1.73m}^2$ or
- FSGS recurs after renal transplantation.

Treatment Frequency:

- **2x** a week for **3** weeks followed by **1x** a week for **6** weeks.
- Each treatment lasts about **2-4** hours.

**Just 12
Treatments
Over
9 Weeks**

How LIPOSORBER LA-15 Works:



1. Blood is withdrawn and goes through the plasma separator.
2. Plasma passes through the LIPOSORBER Adsorption Column, selectively removing LDL, Lp(a) and VLDL.
3. Plasma is recombined with blood cells and returns to the blood warmer, then is returned to the patient.
4. When the primary LIPOSORBER Adsorption Column processes certain amounts of blood, the computer-regulated machine automatically switches the plasma flow to the secondary column.
5. The primary column is regenerated, eluting waste and re-primed to be ready for the next cycle of adsorption.

LIPOSORBER is the only FDA authorized Humanitarian Use Device (HUD) for the treatment of patients with nephrotic syndrome associated with primary FSGS.¹



800-526-3522
liposorber.com


623 Fifth Ave, 27th FL
New York, NY 10022

Kaneka

LIPOSORBER® is a registered trademark of Kaneka Corporation in the United States. When Drug Therapy Fails™ is a claimed common law trademark of Kaneka Corporation in the United States. ©2022 Kaneka Pharma America LLC. All Rights Reserved. BMMKT0002

LIPOSORBER® SELECTIVITY FEATURE

The adsorbent columns are made up of **dextran sulfate cellulose beads**, which provide specific binding to certain lipids such as LDL-C, Lp(a), and VLDL.



**LIPOSORBER
Reduces LDL-C
and Total
Cholesterol³**

By effectively decreasing total cholesterol and LDL-C, LIPOSORBER promotes:

- ➔ **Direct podocytes effects:**
 - Decreases lipotoxicity to the glomeruli⁴
- ➔ **Immunologic effects:**
 - Improves response to corticosteroid and cyclosporine²

Unlike therapeutic plasma exchange (TPE), **LIPOSORBER minimally affects other blood components.**⁵

EFFICACY OF LIPOSORBER




Remission of proteinuria

LIPOSORBER LA-15 therapy has induced complete or partial remission of proteinuria in up to 50% of patients with both drug-resistant and recurrent FSGS.³

Improved responsiveness



Lipoprotein-apheresis (LA) therapy improves corticosteroid and cyclosporine action and responsiveness in patients with drug resistance, which is likely due to the restoration of cellular uptake and the inhibition of drug efflux.²



ADVERSE EVENTS: The most common adverse events are hypotension (0.8%), nausea/vomiting (0.5%), and flushing/blotching (0.4%). Other adverse reactions include angina/chest pain, shortness of breath, fainting, lightheadedness and anemia.¹ *Please see liposorber.com for a full list of adverse events.



CONTRAINDICATION: Angiotensin converting enzyme [ACEI(s)] inhibitors are contraindicated with LIPOSORBER® due to possible bradykinin reaction. ACEI(s) should be replaced with angiotensin II receptor blockers (ARBs) or any other antihypertensive agent as determined by the prescribing physician.¹

References:

1. Kaneka Medical America LLC., 2021, LIPOSORBER®LA-15 SYSTEM Operator's Manual No. 1002en-R4
2. Muso, E. Beneficial effect of LDL-apheresis in refractory nephrotic syndrome. Clin Exp Nephrol 2014; 18: 286-290.
3. Muso, E et al. Immediate therapeutic efficacy of LDL apheresis for drug-resistant nephrotic Syndrome: evidence from the short-term results from the POLARIS Study. Clin Exp Nephrol. 2015; 19: 379-386.
4. Raina R, Wang J, Sharma A, Chakraborty R: Extracorporeal Therapies in the Treatment of Focal Segmental Glomerulosclerosis. Blood Purif 2020;49:513-523. doi: 10.1159/000506277
5. Rubba, Paolo et al. Hemodynamic changes in the peripheral circulation after repeat low density lipoprotein apheresis in familial hypercholesterolemia. Circulation 81 2 (1990): 610-6.